

# Wet Granulation: End-Point Determination

by Michael Levin, Ph. D.

Measurement Control Corporation  
East Hanover, New Jersey, USA

## What is an end-point?

End-point can be defined by the formulator as a target particle size mean or distribution. Alternatively, the end-point can be defined in rheological terms. It has been shown (15) that once you have reached the desired end-point, the granule properties and the subsequent tablet properties are very similar regardless of the granulation processing factors, such as impeller or chopper speed or binder addition rate. I would call this “the principle of equifinality”.

The ultimate goal of any measurement in a granulation process is to estimate viscosity and density of the granules, and, perhaps, to obtain an indication of the particle size mean and distribution. One of the ways to obtain this information is by measuring load on the main impeller.

Mixer instrumentation, in general, has numerous benefits. In addition to a possible end-point determination, it can be used to troubleshoot the machine performance (for example, help detect worn-out gears and pulleys or identify mixing and binder irregularities). Instrumentation can serve as a tool for formulation fingerprinting, assure batch reproducibility, aid in raw material evaluation, process optimization and scale-up.

## End-Point Determination

End-point detection in wet granulation has become a major scientific and technological challenge (74). Monitoring granulation is most commonly achieved by collecting either power or torque signals, or both. In what follows, we will compare both methods.

How to determine an end-point?

A wet granulation end-point should be defined empirically in terms of wet mass density and viscosity, particle size distribution, flowability or tableting parameters (e.g., capping compression).

It is advisable to run a trial batch at a fixed speed and with a predetermined method of binder addition (for example, add water continuously at a fixed rate to a dry mix with a water-soluble binding agent).

Before adding the liquid, measure the baseline level of motor power consumption  $P_o$  or impeller torque  $T_o$  at the dry mix stage.

During the batch, stop the process frequently times to take samples and, for each sample, note the end-point values of power consumption  $P_e$  or impeller torque  $T_e$ . For each of these “end-points”, measure the resulting wet mass density  $\rho$ . As a result, you will be able to obtain some data that will relate the “end-point parameters” listed above with the processing variables in terms of net motor power consumption  $\Delta P_m = (P_e - P_o)$  or net

impeller power consumption  $\Delta P_i = 2\pi \cdot n \cdot (T_e - T_o)$ , where  $n$  is the impeller speed [dimension  $T^{-1}$ ].

## Torque vs. Power

When we say “power consumption”, we usually refer to the main motor. It reflects the load on the motor due to useful work, as well as the power needed to run the motor itself (losses due to eddy currents, friction in couplings, etc.).

It is quite possible (and, indeed, quite pertinent) to talk about the power consumption of the impeller, which is, obviously, quantitatively less than the power consumption of the motor and relates directly to the load on the impeller.

### **Power ~ Torque \* Speed**

Impeller power consumption can be calculated as a product of the direct torque, rotational impeller speed, and a coefficient (usually equal to  $2\pi$  times a unit conversion factor, if required).

The power consumption of the mixer motor differs from that of the impeller by the variable amount of power draw imposed by various sources (mixer condition, transmission, gears, couplings, motor condition, etc.)


Compared to impeller torque, motor power consumption is easier to measure; watt meters are inexpensive and can be installed with almost no downtime. However, motor power signal may not be sensitive enough for specific products or processing conditions. Wear and tear of mixer and motor may cause power fluctuations. Moreover, power baseline may shift with load.

Impeller torque, on the other hand, is closer to where the action is, and is directly related to the load on the impeller. Torque is not affected by mixer condition.

Although the motor power consumption is strongly correlated with the torque on the impeller (38), it is less sensitive to high frequency oscillations caused by direct impact of particles on the blades as evidenced by FFT technique (16).

Power consumption or torque fluctuations are influenced by granule properties (particle size distribution, shape index, apparent density) and the granulation time. Fluctuation of torque / power consumption and intensity of spectrum obtained by FFT analysis can be used for end-point determination (37).

It was observed that when the end-point region of a granulation is reached, the frequency distribution of a power consumption signal reaches a steady state (75). It should be repeated here that torque shows more sensitivity to high frequency oscillations.



### Torque and Power Profiles

Fig. 1. Impeller torque and motor power consumption for a small high shear mixer (Felder PMA 10).

Fig. 1 illustrates the classical power and torque profiles that start with a dry mixing stage, rise steeply with binder solution addition, level off into a plateau, and then exhibit overgranulation stage. The power and torque signals have similar shape and are strongly correlated. The pattern shows a plateau region where power consumption or torque is relatively stable.

The peak of the derivative indicates the inflection point of the signal. Based on the theory by Leuenberger (1979 and subsequent work), useable granulates can be obtained in the region that starts from the peak of the signal derivative with respect to time and extends well into the plateau area (40). Prior to the inflection point, a continuous binder solution addition may require variable quantities of liquid. After that point, the process is well defined and the amount of binder solution required to reach a desired end-point may be more or less constant.

Torque or power consumption pattern of a mixer is a function of the viscosity of both the granulate and binder. With the increasing viscosity, the plateau is shortened and sometimes vanishes completely thereby increasing the need to stop the mixer at the exact end-point.

At low impeller speeds or high liquid addition rates, the classic S-shape of the power consumption curve may become distorted with a steep rise leading into overgranulation (9).

The area under the torque-time curve is related to the energy of mixing and can be used as an end-point parameter. Area under power consumption curve divided by the load gives the specific energy consumed by the granulation process. This quantity is well

correlated with the relative swept volume (11, 12, 32).

The consumed energy is completely converted into heat of the wet mass (7), so that the temperature rise during mixing shows some correlation with relative swept volume and Froude number (29) that relates the inertial stress to the gravitational force per unit area acting on the material.




Fig. 2. A torque profile in a typical production batch

Fig. 2 represents a record of a typical granulation batch done by an experienced operator on large Hobart mixer. You can see that the batch was stopped on the downslope of the derivative.




Fig. 3. Another batch by the same operator (power consumption profile)

On a Fig. 3 you can see another batch made by the same operator. This time it is a power consumption trace, but again it extends beyond the peak of the derivative and the end-point thus can be deemed reproducible.




Fig.4. A batch by a novice operator (power consumption profile)

In the batch represented in Fig. 4, a novice operator trainee has stopped the batch well before the peak of the derivative. This required a major adjustment of the tableting operation (force and speed) to produce tablets in an acceptable range of material properties (hardness and friability).





Fig. 5. Another batch by inexperienced operator (torque profile)

In this batch (Fig. 5), the same novice operator has stopped the granulation process, opened the lid, took a sample, and decided to granulate for another 10 seconds. You can see that there is no indication that the peak of the derivative was reached at the end-point.

Thus, it seems that monitoring torque or power can fingerprint not only the product, but the process and the operators as well.

A number of publications relate to practical experience of operators on the production floor (34, 76-78).



## End-Point Optimization

Fig. 6. Wet granulation end-point as a factor in tableting optimization

Agglomeration of particles in wet granulation have been studied extensively (24, 79). The optimal end-point can be thought of as the factor affecting a number of agglomerate properties (Fig.8).

With so many variables involved in a granulation process, it is no wonder that more and more researchers throw in a number of factors together in an attempt to arrive at an optimum response (30, 35, 80-88).

The final goal of any granulation process is a solid dosage form, such as tablets. Therefore, when optimizing a granulation process, the list of factors affecting tablet properties may include both the granulation end-point and the tableting processing parameters, such as compression force or tablet press speed.

In one of the most interesting works based on the experimental design approach, an attempt was made to find a statistical relationship between the major factors affecting both granulation and compaction, namely, granulation end-point, press speed (dwell time), and compression force (88). The resulting equation allowed optimization of such standard response parameters as tablet hardness, friability and disintegration time. This study has also investigated the possibility of adjusting the tableting parameters in order to account for an inherent variability of a wet granulation process.

Multivariate optimization of wet granulation may include hardness, disintegration and ejection as response variables (89). Compressibility property of granulations is extremely sensitive to various processing parameters of wet granulation (90).



Recently, the experimental design procedure was applied to low shear wet granulation (91) with a factorial design used to evaluate the influence of such factors as binder strength and agitator speed.

## End-Point Reproducibility

As will be shown in the next section, for every blend and a fixed set of values for processing factors (such as mixer geometry, blade speed, powder volume, amount and method of addition of granulating liquid), a wet granulation process state (end-point) is completely characterized by rheological properties of the wet mass (density, viscosity), which are, in turn, a function of particle size, shape and other properties. The process can be quantified with the help of dimensionless Newton Power Number  $N_p$  that will assume a certain numerical value for every state (condition) of the granulate. Under fixed processing conditions,  $N_p$  will be proportional to Net Power Consumption  $\Delta P$  for any end-point (defined, in part, by wet mass density). Thus, in order to reproduce an end-point, it is sometimes sufficient to monitor power of the impeller (or the motor) and stop when a predefined net level of the signal is reached. If, however, any of the processing variables or the rheological definition of the end-point has changed, a more sophisticated approach is required, as described below.

Once the desired end-point is determined, it can be reproduced by stopping the batch at the same level of net power consumption  $\Delta P$  (for the same mixer, formulation, speed, batch size and amount/rate of granulating liquid). To account for changes in any of these variables, you have to compute the Newton power number  $N_p$  for the desired end-point:

$$N_p = \Delta P / (\rho n^3 d^5)$$

In other words, if you have established an end-point in terms of some Net Impeller or Motor Power  $\Delta P$  and would like to reproduce this end-point on the same mixer at a different speed or wet mass density, calculate Newton Power Number  $N_p$  from the given Net Impeller Power  $\Delta P$ , impeller speed  $n$ , blade radius  $d$ , and wet mass density  $\rho$  (assuming the same batch size), and then recalculate the target  $\Delta P$  with the changed values of speed  $n$  or wet mass density  $\rho$ .

Wet mass viscosity  $\eta$  can be calculated from Net Impeller Power  $\Delta P$ , blade radius  $d$  and impeller speed  $n$ , using the following equations:

$$\begin{aligned} \Delta P &= 2\pi \Delta\tau * n \\ \eta &= \phi * \Delta\tau / (n * d^3) \end{aligned}$$

where  $\Delta\tau$  is the net torque required to move wet mass,  $n$  is the speed of the impeller,  $d$  is the blade radius or diameter, and  $\phi$  is mixer specific “viscosity factor” relating torque and dynamic viscosity (note: the correlation coefficient  $\phi$  can be established empirically by mixing a material with a known dynamic viscosity, e.g. water). Alternatively, you can use impeller torque  $\tau$  as a measure of kinematic viscosity and use it to obtain a non-dimensionless “pseudo-Reynolds” number, based on the so-called “mix consistency” measure, that is, the end-point torque, as described in the case studies.

Fill Ratio  $h/d$  can be calculated from a powder weight, granulating liquid density (1000 kg/m<sup>3</sup> for water), rate of liquid addition, time interval for liquid addition, and bowl volume  $V_b$ . The calculations are performed using the idea that the fill ratio  $h/d$  (wet mass height to blade diameter) is proportional to  $V/V_b$ , and wet mass volume  $V$  can be computed as

$$V = m / \rho,$$

where  $m$  is the mass (weight) of the wet mass and  $\rho$  is the wet mass density.

Now, the weight of the wet mass is computed as the weight of powder plus the weight of added granulating liquid. The latter, of course, is calculated from the rate and duration of liquid addition and the liquid density.

Finally, you can combine the results obtained at different end-points of the test batch or from different batches or mixer scales (assuming geometrical similarity).

Given wet mass density  $\rho$ , wet mass viscosity  $\eta$ , fill ratio  $h/d \sim m/V_b / \rho$ , setup speed  $n$ , and blade radius or diameter  $d$ , you can calculate the Reynolds number  $Re$  (or the "pseudo-Reynolds" number) and the Froude number  $Fr$ . Then you can estimate the slope "a" and intercept "b" of the regression equation

$$N_p = b \cdot (Re \cdot Fr \cdot h/d)^a$$

or

$$\log N_p = \log b + a \cdot \log (Re \cdot Fr \cdot h/d)$$

And, inversely, once the regression line is established, you can calculate Newton Power number  $N_p$  (which is the target quantity for scale up) and Net Power  $\Delta P$  (which can be observed in real time as a true indicator of the target end-point) for any point on the line (40, 42, 100-103, ).

# List of Symbols and Dimensions

a, b	Slope and intercept of a regression equation
d	Impeller (blade) diameter or radius (m); dimensional units [L]
g	Gravitational constant (m / s <sup>2</sup> ); dimensional units [LT <sup>-2</sup> ]
h	Height of granulation bed in the bowl (m); dimensional units [L]
H	Bowl height (m); dimensional units [L]
l	Blade length (m); dimensional units [L]
n	Impeller speed (revolutions / s); dimensional units [T <sup>-1</sup> ]
P	Power required by the impeller or motor (W = J / s); dimensional units [ML <sup>2</sup> T <sup>-5</sup> ]
R <sub>b</sub>	Radius of the bowl (m); dimensional units [L]
q	Binder liquid addition rate
s	Amount of granulating liquid added per unit time (kg); dimensional units [M]
t	Binder addition time (s); dimensional units [T]
V <sub>p</sub>	Particle volume (m <sup>3</sup> ); dimensional units [L <sup>3</sup> ]
V <sub>m</sub>	Wet mass volume (m <sup>3</sup> ); dimensional units [L <sup>3</sup> ]
V <sub>b</sub>	Bowl volume (m <sup>3</sup> ); dimensional units [L <sup>3</sup> ]
w	Wet mass; dimensional units [M]
ρ	Specific density of particles (kg / m <sup>3</sup> ); dimensional units [M L <sup>-3</sup> ]
$\nu = \eta / \rho$	Kinematic viscosity (m <sup>2</sup> / s); dimensional units [L <sup>2</sup> T <sup>-1</sup> ]
η	Dynamic viscosity (Pa*s); dimensional units [M L <sup>-1</sup> T <sup>-1</sup> ]
τ	Torque (N-m); dimensional units [M L <sup>2</sup> T <sup>-2</sup> ]. End point torque values were described as “wet mass consistency” numbers. Note: torque has the same dimensions as work or energy.
$\phi = \eta \cdot n \cdot d^3 / \Delta\tau$	Dimensionless “viscosity factor” relating net torque Δτ and dynamic viscosity η
$Fr = n^2 d / g$	Froude number. It relates the inertial stress to the gravitational force per unit area acting on the material. It is a ratio of the centrifugal force to the gravitational force.
$N_p = P / (\rho n^3 d^5)$	Newton (power) number. It relates the drag force acting on a unit area of the impeller and the inertial stress.
$Re = d^2 n \rho / \eta$	Reynolds number. It relates the inertial force to the viscous force.
$\Psi Re = d^2 n \rho / \tau$	“Pseudo Reynolds number” (m <sup>3</sup> /s); dimensional units [L <sup>-3</sup> T]. Note: this variable physically is a reciprocal of volume flow rate.
$Ga = Re^2 / Fr$	Galileo number

# Literature Reference

## Reference

- 1 Parikh D. Handbook of Pharmaceutical Granulation Technology, Marcel Dekker, Inc. New York, 1997.
- 2 Sheskey PJ, Williams DM. Comparison of low-shear and high-shear wet granulation techniques and the Influence of percent water addition in the preparation of a controlled-release matrix tablet containing HPMC and a high-dose, highly water-soluble drug. Pharm Tech 3:80-92, 1996
- 3 Morris KR, Nail SL, Peck GE, Byrn SR, Griesser UJ, Stowell JG, Hwang SJ, Park K. Advances in pharmaceutical materials and processing. Pharm Sci Technol Today 1(6):235-245, 1998
- 4 Hausman DS. Comparison of Low Shear, High Shear, and Fluid Bed Granulation During Low Dose Tablet Process Development. Drug Dev Ind Pharm, 30(3):259-266, 2004.
- 5 Holm P, Jungersen O, Schæfer T, Kristensen HG. Granulation in high speed mixers. Part I: Effect of process variables during kneading. Pharm Ind 45:806-811, 1983
- 6 Holm P, Jungersen O, Schæfer T, Kristensen HG. Granulation in high speed mixers. Part II: Effect of process variables during kneading. Pharm Ind 46:97-101, 1984
- 7 Holm P, Schaefer T, Kristensen HG. Granulation in high speed mixers. Part V: Power consumption and temperature changes during granulation. Powder Technol 43:213-223, 1985
- 8 Holm P, Schaefer T, Kristensen HG. Granulation in high-speed mixers. Part IV. Effects of process conditions on power consumption and granule growth. Powder Technol 43:225, 1985
- 9 Holm P, Schaefer T, Kristensen HG. Granulation in high speed mixers. Part VI: Effects of process conditions on power consumption and granule growth. Powder Technol 3:286, 1993
- 10 Jaegerskou A, Holm P, Schæfer T and Kristensen HG. Granulation in high speed mixers. Part III: effects of process variables on intergranular porosity. Pharm Ind 46:310-314, 1984
- 11 Schaefer T, Bak HH, Jaegerskou A, Kristensen A, Svensson JR, Holm P and Kristensen HG. Granulation in different types of high speed mixers. Part 1: Effects of process variables and up-scaling. Pharm Ind 48:1083, 1986
- 12 Cliff MJ. Granulation end-point and automated process control of mixer-granulators: Part 1. Pharm Tech 4:112-132, 1990
- 12 Schaefer T, Bak HH, Jaegerskou A, Kristensen A, Svensson JR, Holm P and Kristensen HG. Granulation in different types of high speed mixers. Part 2: Comparison between mixers. Pharm Ind 49:297-304, 1987
- 13 Cliff MJ. Granulation end-point and automated process control of mixer-granulators: Part 2. Pharm Tech 5:38-44, 1990

- 15 Emori H, Sakuraba Y, Takahashi K, Nishihata T, Mayumi T. Prospective validation of high-shear wet granulation process by wet granule sieving method. II. Utility of wet granule sieving method. *Drug Dev Ind Pharm*, 23(2):203-215, 1997
- 16 Corvari V, Fry W C, Seibert WL, Augsburger L. Instrumentation of a high-shear mixer: Evaluation and comparison of a new capacitive sensor, a watt meter, and a strain-gage torque sensor for wet granulation. *Pharm Res* 9(12):1525-1533, 1992
- 17 Corvari V, Fry W C, Seibert WL, Augsburger L. Wet granulation end-point detection in a high shear mixer instrumented with a capacitive sensor and a strain gaged torque sensor. AAPS Meeting, 1992
- 18 Fry WC, Stagner WC, Wichman KC. Computer-interfaced capacitive sensor for monitoring the granulation process 1: Granulation monitor design and application. *J Pharm Sci* 73:420-421, 1984
- 19 Fry WC, Stagner WC, Wichman KC. Computer-interfaced capacitive sensor for monitoring the granulation process 2: System response to process variables. *Pharm Tech* 30-41, Oct, 1987
- 20 Terashita K, Kato M, Ohike A. Analysis of end-point with power consumption in high speed mixer. *Chem Pharm Bull* 38(7):1977-1982, 1990
- 21 Spring MS. The conductivity of the damp mass during the massing stage of the granulation process. *Drug Dev Ind Pharm* 9(8), 1507-1512, 1983
- 22 Staniforth JN, Quincey SM. Granulation monitoring in a planetary mixer using a probe vibration analysis technique. *Int J Pharm* 32, 177-185, 1986
- 23 Kay D, Record PC. Automatic wet granulation end-point control system. *Manuf Chem Aerosol News* 9:45-46, 1978
- 23 Staniforth JN, Walker S, Flander P. Granulation monitoring in a high speed mixer/processor using a probe vibration analysis technique. *Int J Pharm* 31, 277-280, 1986
- 24 Alderborn G. Granule properties of importance to tableting. *Acta Pharm Scand* 25:229-238, 1988
- 25 Timko RJ, Johnson JL, Skinner GW, Chen ST, Rosenberg HA. Instrumentation of a vertical high shear mixer with a motor slip monitoring device. *Drug Dev Ind Pharm* 12(10):1375-1393, 1986
- 26 Timko RJ, Barrett JS, McHugh PA, Chen ST, Rosenberg HA. Use of a motor load analyzer to monitor the granulation process in a high intensity mixers. *Drug Dev Ind Pharm* 13(3):405-435, 1987
- 27 Fink DG, Beaty HW. *Standard Handbook for Electrical Engineers*. 13th edition. McGraw-Hill, New York 2-17, 3-26:27, 20-13, 20-40, 1993
- 28 Oldshue JY. Mixing processes. In: Bisio A, Kabel RL (eds). *Scale-up of Chemical Processes: Conversion from Laboratory Scale Tests to Successful Commercial Size Design*. Wiley, New York, 1985

- 29 Horsthuis GJB, van Laarhoven JAH, van Rooij RCBM, Vromans H. Studies on upscaling parameters of the Gral high shear granulation process. *Int J Pharm* 92:143, 1993
- 30 Rekhi GS, Caricofe RB, Parikh DM, Augsburger L L. A new approach to scale-up of a high-shear granulation process. *Pharm Tech Suppl - TabGran Yearbook* 58-67, 1996
- 31 Schaefer T. Equipment for wet granulation. *Acta Pharm Seuc*, 25:205, 1988
- 32 Holm P. Effect of impeller and chopper design on granulation in a high speed mixer. *Drug Dev Ind Pharm* 13:1675, 1987
- 34 Werani J. Production experience with end-point control. *Acta Pharm Seuc* 25:247-266, 1988
- 35 Lindberg N-O, Jonsson C, Holmquist B. The granulation of a tablet formulation in a high-speed mixer, Diosna P25. *Drug Dev Ind Pharm* 11:917-930, 1985
- 36 Laicher A, Profitlich T, Schwitzer K, Ahlert D. A modified signal analysis system for end-point control during granulation. *Eur J Pharm Sci* 5:7-14, 1997
- 37 Kristensen HG, Holm P, Jaegerskou A and Schaefer T. Granulation in high speed mixers. Part IV: Effect of liquid saturation on the agglomeration. *Pharm Ind* 46:763-767, 1984
- 37 Watano S, Terashita K, Miyanami K. Frequency analysis of power consumption in agitation granulation of powder materials with sparingly soluble acetaminophen. *Chem Pharm Bull* 40(1):269-271, 1992
- 38 Kopcha M, Roland E, Bubb G, Vadino WA. Monitoring the granulation process in a high shear mixer/granulator: an evaluation of three approaches to instrumentation. *Drug Dev Ind Pharm* 18(18):1945-1968, 1992
- 38 Yliruusi J, Tihtonen R. High-speed granulation. Part I: Effect of some process parameters of granulation on the properties of unsieved and wet-sieved granules. *Acta Pharm Fen* 98:39, 1989
- 39 Stamm A, Paris L. Influence of technological factors on the optimal granulation liquid requirement measured by power consumption. *Drug Dev Ind Pharm* 11(2&3), 330-360, 1985
- 39 Yliruusi J, Tihtonen R. High-speed granulation. Part II: Effect of some process parameters on the properties of wet- and dry-sieved granules. *Acta Pharm Fen* 98:53, 1989
- 40 Leuenberger H, Bier HP, Sucker HB. Theory of the granulating-liquid requirement in the conventional granulation process. *Pharm Tech* 6:61-68, 1979
- 40 Schwartz JB, Szymczak CE. Power consumption measurements and the mechanism of granule growth in a wet granulation study. *AAPS Meeting* November, 1997
- 41 Landin M, Rowe RC, York P. Characterization of wet powder masses with a mixer torque rheometer. 3. Nonlinear effects of shaft speed and sample weight. *J Pharm Sci* 84/5:557-560, 1995

- 42 Landin M, York P, Cliff MJ, Rowe RC, Wigmore AJ. The effect of batch size on scale-up of pharmaceutical granulation in a fixed bowl mixer-granulator. *Int J Pharm* 134:243-246, 1996
- 43 Landin M, York P, Cliff MJ. Scale-up of a pharmaceutical granulation in fixed bowl mixer granulators. *Int J Pharm* 133:127-131, 1996
- 44 Landin M, York P, Cliff MJ, Rowe RC. Scaleup of a pharmaceutical granulation in planetary mixers. *Pharm Dev Technol*, 4(2):145-150, 1999
- 45 Faure A, Grimsey IM, Rowe RC, York P, Cliff MJ. A methodology for the optimization of wet granulation in a model planetary mixer. *Pharm Dev Tech* 3(3):413-422, 1998
- 46 Faure A, Grimsey IM, Rowe RC, York P, Cliff MJ. Importance of wet mass consistency in the control of wet granulation by mechanical agitation: a demonstration. *J Pharm Pharmacol*, 50(12):1431-2, 1998
- 47 Faure A, Grimsey IM, Rowe RC, York P, Cliff MJ. Applicability of a scale-up methodology for wet granulation processes in Collette Gral high shear mixer-granulators. *Eur J Pharm Sci*, 8(2):85-93, 1999
- 48 Faure A, Grimsey IM, Rowe RC, York P, Cliff MJ. Process control in a high shear mixer-granulator using wet mass consistency: The effect of formulation variables. *J Pharm Sci*, 88(2):191-195, 1999
- 49 Faure A, York P, Rowe RC. Process control and scale-up of pharmaceutical wet granulation processes: a review. *Eur J Pharm Biopharm*, 52(3):269-277, 2001
- 50 Betz G, Bürgin PJ, Leuenberger H. Power consumption profile analysis and tensile strength measurements during moist agglomeration. *Int J Pharm*, 252(1-2):11-25, 2003
- 51 Betz G, Bürgin PJ, Leuenberger H. Power consumption measurement and temperature recording during granulation. *Int J Pharm*, 272(1-2):137-149, 2004
- 52 Ritala M, Holm P, Schaefer T, Kristensen HG. Influence of liquid bonding strength on power consumption during granulation in a high shear mixer. *Drug Dev Ind Pharm* 14:1041, 1988
- 53 Sirois PJ, Craig GD. Scaleup of a high-shear granulation process using a normalized impeller work parameter. *Pharm Dev Technol* 5(3):365-374, 2000
- 54 Hirzel J. Understanding premium-efficiency motor economics. *Plant Eng* May 7:75-78, 1992
- 55 Elliott T. Efficiency, reliability of drive systems continue to improve. *Power* February:33-41, 1993
- 56 Holm P. High shear mixer granulators. In: Parikh DM (ed.). *Handbook of Pharmaceutical Granulation Technology*. Marcel Dekker, Inc. New York, 1997
- 57 Chirkot T, Propst CW. Low shear granulators. In: Parikh DM (ed.). *Handbook of Pharmaceutical Granulation Technology*. Marcel Dekker, Inc. New York, 1997
- 58 Ghanta SR, Srinivas R, Rhodes CT. Use of mixer-torque measurements as an aid to optimizing wet granulation process. *Drug Dev Ind Pharm* 10(2):305-311, 1984



- 59 Rowe RC, Sadeghnejad GR. The rheological properties of microcrystalline cellulose powder/water mixes - measurement using a mixer torque rheometer. *Int J Pharm* 38:227-229, 1987
- 60 Hancock BC, York P, Rowe RC. Characterization of wet masses using a mixer torque rheometer: 1: Effect of instrument geometry. *Int J Pharm* 76:239-245, 1991
- 61 Hariharan M, Mehdizadeh E. The Use of Mixer Torque Rheometry to Study the Effect of Formulation Variables on the Properties of Wet Granulations. *Drug Dev Ind Pharm*, 28(3):253-263, 2002
- 62 Parker MD, Rowe RC, Upjohn NG. Mixer torque rheometry: A method for quantifying the consistency of wet granulation's. *Pharm Tech Int* 2:50-64, 1990
- 63 Rowe RC, Parker MD. Mixer torque rheometry: An update. *Pharm Tech* 74-82, March, 1994
- 64 Watano S, Sato Y, Miyanami K. Application of a neural network to granulation scale-up. *Powder Technol*, 90(2):153-159, 1997
- 65 Watano S. Direct control of wet granulation processes by image processing system. *Powder Technol*, 117(1-2):163-172, 2001
- 66 Laurent BFC. Structure of powder flow in a planetary mixer during wet-mass granulation. *Chem Eng Sci*, 60(14):3805-3816, 2005
- 67 Whitaker M, Baker GR, Westrup J, Goulding PA, Rudd DR, Belchamber RM, Collins MP. Application of acoustic emission to the monitoring and end-point determination of a high shear granulation process. *Int J Pharm*, 205(1-2):79-92, 2000
- 68 Belchamber R. Acoustics - a process analytical tool. *Spec Eur*, 15(6):26-7, 2003
- 69 Rudd D. The Use of Acoustic Monitoring for the Control and Scale-Up of a Tablet Granulation Process. *J Proc Anal Tech*, 1(2):8-11, 2004
- 70 Miwa A, Toshio Yajima T, Itai S. Prediction of suitable amount of water addition for wet granulation. *Int J Pharm*, 195(1-2):81-92, 2000
- 71 Otsuka M, Mouri Y, Matsuda Y. Chemometric Evaluation of Pharmaceutical Properties of Antipyrine Granules by Near-Infrared Spectroscopy. *AAPS Pharm Sci Tech*, 4(3) Article 47, 2003
- 72 Ganguly S, Gao JZ. Application of On-line Focused Beam Reflectance Measurement Technology in High Shear Wet Granulation. *AAPS General Meeting, Contributed Paper*, 2005
- 73 Dilworth SE, Mackin LA, Weir S, Claybourn M, Stott PW. In-line techniques for end-point determination in large scale high shear wet granulation. *142nd British Pharmaceutical Conference*, 2005.
- 74 Holm P, Schaefer T, Larsen C. End-Point Detection in a Wet Granulation Process. *Pharm Dev Technol*, 6(2):181-192, 2001
- 75 Terashita K, Watano S, Miyanami K. Determination of end-point by frequency analysis of power consumption in agitation granulation. *Chem Pharm Bull* 38(11):3120-3123, 1990



- 76 Titley PC. Agglomeration and granulation of powders, processing and manufacturing practice. *Acta Pharm Seuc* 25:267-280, 1988
- 77 Lindberg N-O. Some experience of continuous granulation. *Acta Pharm Seuc*. 25:239-246, 1988
- 78 Record PC. Practical experience with high-speed pharmaceutical mixer/granulators. *Manuf Chem Aerosol News* 11:65, 1979
- 79 Kristensen HG. Agglomeration of powders. *Acta Pharm Seuc* 25:187-204, 1988
- 80 Vojnovic D, Selenati P, Rubessa F, Moneghini M. Wet granulation in a small scale high shear mixer. *Drug Dev Ind Pharm* 18:961, 1992
- 81 Wehrle P, Nobelis P, Cui n  A, Stamm A. Response surface methodology: an interesting statistical tool for process optimization and validation: example of wet granulation in a high-shear mixer. *Drug Dev Ind Pharm* 19:1637, 1993
- 82 Vojnovic D, Moneghini M, Rubessa F. Simultaneous optimization of several response variables in a granulation process. *Drug Dev Ind Pharm* 19:1479, 1993
- 83 Vojnovic D, Moneghini M, Rubessa F. Optimization of granulates in a high shear mixer by mixture design. *Drug Dev Ind Pharm* 20:1035, 1994
- 84 Miyamoto Y, Ogawa S, Miyajima M, Sato H, Takayama K, Nagai T. An evaluation of process variables in wet granulation. *Drug Dev Ind Pharm* 21:2213, 1995
- 85 Miyamoto Y, Ogawa S, Miyajima M, Matsui M, Sato H, Takayama K, Nagai T. An application of the computer optimization technique to wet granulation process involving explosive growth of particles. *Int J Pharm*, 149(1):25-36, 1997
- 86 Miyamoto Y, Ryu A, Sugawara S, Miyajima M, Ogawa S, Matsui M, Takayama K, Nagai T. Simultaneous optimization of wet granulation process involving factor of drug content dependency on granule size. *Drug Dev Ind Pharm*, 24(11):1055-1056, 1998
- 87 Ogawa S, Kamijima T, Miyamoto Y, Miyajima M, Sato H, Takayama K, Nagai T. A new attempt to solve the scale-up problem for granulation using response surface methodology. *J Pharm Sci* 83(3):439-443, 1994
- 88 Iskandarani B, Shiromani PK, Clair JH. Scale-up Feasibility in High-Shear Mixers: Determination Through Statistical Procedures. *Drug Dev Ind Pharm*. 27(7):651-657, 2001
- 89 Achanta AS, Adusumilli P, James KW. end-point determination and its relevance to physicochemical characteristics of solid dosage forms. *Drug Dev Ind Pharm* 23(6):539-546, 1997
- 89 Westerhuis JA, Coenegracht PMJ, Lerk CF. Multivariate modelling of the tablet manufacturing process with wet granulation for tablet optimization and in-process control. *Int J Pharm* 156(1):109-117, 1997
- 90 Badawy SIF, Menning M, Gorko MA, Gilbert DL. Effect of process parameters on compressibility of granulation manufactured in a high-shear mixer. *Int J Pharm* 198(1):51-61, 2000
- 91 Chirkot T. Scale-Up and Endpoint Issues of Pharmaceutical Wet Granulation in a V-Type Low Shear Granulator. *Drug Dev Ind Pharm*. 28(7):871-888, 2002

- 93 Zega J, Lee D, Shiloach A, Erb D. Scale-up of the wet granulation process for a dicalcium phosphate formulation using impeller power consumption. AAPS Meeting November, 1995
- 94 Rayleigh Lord. The principle of similitude. *Nature* 95(2368, March 18):66-68, 1915
- 95 Leuenberger H. Granulation, new technique. *Pharm Acta Helv* 57(3):72-80, 1982
- 96 Merrifield CW. The experiments recently proposed on the resistance of ships. *Trans Inst Naval Arch (London)* 11:80-93, 1870
- 96 Zlokarnik M. Problems in the application of dimensional analysis and scale-up of mixing operations. *Chem Eng Sci* 53(17):3023-3030, 1998
- 97 Reynolds O. An experimental investigation of the circumstances which determine whether the motion of water shall be direct or sinusous, and of the law of resistance in parallel channels. *Philos Trans R Soc London* 174:935-982, 1883
- 98 Zlokarnik M. Dimensional analysis and scale-up in chemical engineering. Springer Verlag, 1991
- 99 Buckingham E. On physically similar systems; Illustrations of the use of dimensional equations. *Phys Rev NY* 4:345-376, 1914
- 100 Bier HP, Leuenberger H, Sucker H. Determination of the uncritical quantity of granulating liquid by power measurements on planetary mixers. *Pharm Ind* 4:375-380, 1979
- 101 Leuenberger H. Monitoring granulation. *Manuf Chem Aerosol News* 67-71, May, 1983
- 102 Leuenberger H. Monitoring granulation, Part 2. *Manuf Chem Aerosol News*, June, 1983
- 103 Leuenberger H. Scale-up of granulation processes with reference to process monitoring. *Acta Pharm Technol* 29(4), 274-280, 1983
- 104 Hutin S, Chamayou A, Avan JL, Paillard B, Baron M, Couarraze G, Bougaret J. Analysis of a Kneading Process to Evaluate Drug Substance–Cyclodextrin Complexation. *Pharm Tech.*, pp. 112-123, October 2004