



**Measurement Control Corporation**  
*Process Analytical Instrumentation, Monitoring and Control*  
13D Great Meadow Lane, East Hanover, New Jersey 07936, U.S.A.  
Tel: (800) 504-9010 Fax: (862) 250-6299 email: email@mcc-online.com  
Web Address: [www.mcc-online.com](http://www.mcc-online.com)

# Compression Force Meter for all your Production Tablet Presses!

Compression Force Meter (CFM) from MCC

Inexpensive device designed to monitor force on any tablet press, with high-speed accurate real time display of peak force in kg, kN, or pounds.

It is ideal for press set-up and monitoring, to prevent excessive force.

## The CFM includes:

- Calibrated force transducer
- Maximum force readout
- Signal conditioning
- User setup alarms
- Installation instructions
- Operation instructions
- One year warranty



The meter meets NEMA 4X/IP65 requirements.

**CFM** will display maximum force, with manual pushbutton reset. User can set a force level that will trigger an audible alarm.